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- A system that allows the easy and safe luxation of retained roots that have suffered great damage or are broken off – without surgical intervention
- There is no soft or hard tissue damage affecting the gums and alveolar process surrounding the retained root
- Safe handling: unlike in the case of traditional root luxations that are carried out using elevators, pliers, scalpels or burs, there is no risk of injuring the patient or oneself
- Easy removal of dental germs from bone lacunae that are difficult to inspect and hard to access (germectomy)
- Relaxed ‘extraction experience’, which is not overly stressful for the patient
- Option to carry out immediate implantation after the removal of the retained roots is not impaired
- Due to the special design, this is a safe and easy-to-handle tool that can be taken apart and cleaned quickly

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¹ 205=RA L ² Largest working part diameter in 1/10 mm Teflon® is a registered trademark of DuPont, USA

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
General instructions:

- All products delivered are unsterile therefore, before initial and each further application to be treated (cleaning / disinfection / sterilisation).

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



¹ 204=RA ² Largest working part diameter in 1/10 mm Teflon® is a registered trademark of DuPont, USA

	Wurzelanker Root anchors Ancla de raíz		Extraktionsanker Extraction anchors Ancla de extracción illustrated 1:2		Ratsche Ratchet Carraca illustrated 1:3
					
Fig.	AP23	AP34	AP23	AP34	CARA2
Size ²	016	021	016	021	
Length mm	8,0	8,0	12,0	12,0	

² Largest working part diameter in 1/10 mm

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

85FL072-0609

Hager & Meisinger GmbH

Hansemannstr. 10
41468 Neuss • Germany
Tel.: +49 (0) 21 31-20 12 0
Fax: +49 (0) 21 31-20 12 222
Internet: www.meisinger.de
www.bone-management.de
E-mail: info@meisinger.de

Meisinger USA, L.L.C.

7442 South Tucson Way, Suite 130
Centennial, Colorado 80112 • USA
Tel.: +1 (303) 268-5400
Toll free: +1 (866) 634-7464
Fax: +1 (303) 268-5407
Internet: www.meisingerusa.com
www.bone-management.com
E-mail: info@meisingerusa.com

MEISINGER
GERMANY / USA

ApEx-Control developed with Dr. Guse, MSc

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With the help of ApEx-Control and ApEx-Control Professional, very damaged, retained roots are luxated easily, gently and safely. Compared with the conventional surgical removal of retained roots where gums are pushed aside and milling weakens some of the precious crestal alveolar bone, there is hardly any soft and hard tissue damage when using ApEx-Control. Accidental injuries to the patient or oneself, which, occasionally, may even happen to an experienced surgeon, are a thing of the past. It is unlikely that the patient will experience any pain or swellings. The wound heals much more quickly. Overall, ApEx-Control helps to preserve and maintain desmodontal structures, which, particularly in the case of immediate implantations – which are the preferred procedure - need to be as intact as possible. With the help of ApEx-Control Professional even dental germs that are of a more round shape and therefore difficult to get hold of with the elevator and pliers, can be removed from a bone cavity that is difficult to inspect and hard to access. Furthermore, the Bone Management® System ApEx-Control Professional includes extraction screws of different sizes, shorter pilot drills, a larger spacer and additional extraction anchors.

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ApEx-Control / ApEx-Control Professional

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Clear the surgical site of any soft and hard fragments located on or around the retained root. Cut the fibre system close to the surface (desmodontium) by using the delicate blade of a periotome or desmotome. In the same way, radicular retentions with multiple roots should initially be separated with a bur.

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Removal of retained roots using a support plate

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First, use moderate pressure and water-cooling to drill a pilot hole parallel to the longitudinal axis of the root until you reach the end of the pilot drill blade.

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
Slowly insert a suitable root anchor with attached insertion device and ratchet until the thread disappears completely into the pilot borehole. Afterwards, remove ratchet and insertion device.

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
Align the support plate with the spacers of various heights with the neighbouring teeth in such a way that the support plate runs parallel to the occlusion level with the extraction anchor protruding centrally through its opening. If the retained root that is to be removed is not in line with the neighbouring teeth i.e. if you are extracting a canine tooth, then the spacers have to be removed and placed in an angle.

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
If the retained root is terminal, the support plate is placed on the neighbouring tooth with the short spacer and on the mandibular ridge with the tall spacer.

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Next, insert the bore-aid for root anchor, which, in some cases, may need to be covered with a distance sleeve of corresponding length, through the opening of the support plate and screw manually onto the root anchor.

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Now place the ratchet onto the boring tool. Use one hand to prevent the support plate from slipping and use the other to slowly rotate the ratchet clockwise. In the case of trapped roots, pause for 5 to 10 seconds after having turned the ratchet by around 30 to 60 degrees in order to give the Scharpey´s fibres of the periodontal ligament time to tear. Only then continue rotating. The luxation finally ends as soon as the resistance that you encounter after the initial turning suddenly recedes and the entire retained root appears before you.

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Removal of retained roots without using a support plate

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
Slowly insert a suitable extraction anchor with the ratchet attached until the thread disappears into the pilot borehole as completely as possible. Afterwards, remove the ratchet. A deep-seated multi-rooted tooth should be divided before extraction with ApEx-Control.

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The remaining root is carefully luxated with the type of pliers that are customarily used to loosen and take hold of retained roots. If the thread of the extraction anchor were to pull out, use an extraction anchor that is one size larger. If possible do not enlarge the pilot borehole any further. If the root split lengthwise, it is easy to remove any fragments subsequently.

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