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Crest-Control and Crest-Control Professional contain:

Crest-Control y Crest-Control Professional contienen:

	Diamantierte Kugel Diamond-headed bur Bola diamantada	Osteotomieinstrumente Osteotomy instruments Fresas para osteotomia	Vorkörner Initial bur Fresa inicial	Diamantierte Sägen Diamond saws Sierras diamantadas	Adapter Adaptors Adaptadores	Ratsche Ratchet Carraca					
Fig.	802L	HM331L	859L	859	186RF	231DC	231DC	231DC	SW0KI	SW0LI	CARA1
Shank ¹	205	205	204	204	204	204	204	204			
Size ²	014	010	010	018	018	070	100	130			
Length mm	8,0	5,5	12,0	10,0	12,0	0,3	0,3	0,3			

¹ 204=RA, 205=RA L, ² Largest working part diameter in 1/10 mm

Horizontal spreaders in Crest-Control

Horizontal-Spreader
Horizontal spreaders
Spreader horizontal
(3 Stck. / 3 p. / 3 pzs.)

Fig.	VCD50
Width	6,0
Length	9,75
Height	10,6

Horizontal spreaders in Crest-Control Professional

Horizontal-Spreader Horizontal spreaders
Spreader horizontal
(je 2 Stck. / 2 p. of each size / 2 pzs. c/u)


Fig.	VCD35	VCD50	VCD85
Width	3,5	6,0	8,5
Length	9,75	9,75	9,75
Height	10,6	10,6	10,6

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- Controlled spreading of horizontally resorbed bone by Bone Spreading
- Effective in increasing the width of narrow alveolar ridges
- Controlled, minimally invasive treatment concept
- Suitable for preliminary implant cavity preparation of all current dental implant systems

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Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

85FL039-0609

Hager & Meisinger GmbH

Hansemannstr. 10
41468 Neuss • Germany
Tel.: +49 (0) 21 31-20 12 0
Fax: +49 (0) 21 31-20 12 222
Internet: www.meisinger.de
www.bone-management.de
E-mail: info@meisinger.de

Meisinger USA, L.L.C.

7442 South Tucson Way, Suite 130
Centennial, Colorado 80112 • USA
Tel.: +1 (303) 268-5400
Toll free: +1 (866) 634-7464
Fax: +1 (303) 268-5407
Internet: www.meisingerusa.com
www.bone-management.com
E-mail: info@meisingerusa.com

MEISINGER
GERMANY / USA

Crest-Control Crest-Control Professional

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Crest-Control is a Horizontal Bone Spreading System which allows for the pre-determined and minimal invasive horizontal widening of the alveolar ridge in particular in the distal lower jaw. With the completely new horizontal spreader the alveolar ridge can be expanded up to a maximum of 5 mm. Subsequently, an implant of all current systems can still be used during the mineralisation phase in the widened segment from which the osteointegration due to a highly-potent body growth factor of the existing callus may be accelerated. The Crest-Control Professional contains two additional Horizontal-Spreaders in widths of 3.5 mm and 8.5 mm.

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
Indications for the use of the Crest-Control and Crest-Control Professional are horizontally resorbed sections e.g. the lower jaw - alveolar ridge, which do not allow the placement of a conventional implant without suitable Bone Management®. In addition, Crest-Control is especially suited when placing several implants side by side.

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
A crestal scoring on the alveolar ridge is carried out first using the diamond spherical instrument. This prevents the following diamond saw from slipping.

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
The diamond saw allows a crestal osteotomy to be created along the scored ridge. This enables the diamond instrument to be used to deepen the crestal osteotomy.

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
In the area of the implant opening, with the help of the osteotomy instrument, crestal cut can be deepened and prepared accordingly for the application of the horizontal spreader.

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
Following this the horizontal spreader is positioned in the crestal cut.

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
With the suitable screwdriver the horizontal spreader is slowly opened thereby widening the osteotomy site. In doing so the bone is at the same time horizontally condensed.

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

The now widened implant site enables the insertion of the suitable implant. In addition, the bone strength clearly due to the condensation offers an optimized primary stability.

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A further application possibility to the Horizontal Augmentation with the Crest-Control is found in the implementation according to the biological protocol of the Callus-Distraktion.

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