Guided-Drill-Control

CE 0044

	Bohrhülsen Drill sleeve Casquillo portabrocas	Tiefenbohrer Depth-burs Fresa de profundidad
Fig.	GDH	GDH
Pieces	10	5
Shank ¹		204
Size ²		030
Length mm	6,0	10,0

¹ 204=RA, 205=RA L

² Largest working part diameter in 1/10 mm

Art.-No.: BGD00

	Vorkörner Initial bur Fresa inicial	Allportbohrer Allport burs Fresa „Allport“	Pilotbohrer Pilot burs Fresas piloto
Fig.	186RF	141RF	141RF
Shank ¹	204	204	204
Size ²	018	023	032
Length mm	12,0		15,0

¹ 204=RA, 205=RA L, 206=RA XL

² Largest working part diameter in 1/10 mm

Art.-No.: BGL00

	Führungsbohrer Guide drills Fresa-guía	Tiefenbohrer Depth burs Fresa de profundidad
Fig.	GD2	GD3
Shank ¹	204	204
Size ²	035	035
Length mm	15,0	15,0

¹ 204=RA, 205=RA L, 206=RA XL

² Largest working part diameter in 1/10 mm

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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
Guided-Drill-Control is the easiest alternative to traditional tube-in-tube systems. After producing 3D images (DVT), the dental laboratory designs and produces a drilling template to ensure the safety of the upcoming pilot drilling. Inserting the drill sleeve is an easy and fast process due to the matching sleeve drills (Guided-Drill-Control Lab).

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
As a first step, the pilot drilling is carried out with the pilot drill that perfectly matches to the drill sleeves. During this process, the drive mechanism limits and secures the depth of the bore.

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
Afterwards, depending on the diameter of the planned implant and the quality of the jaw bone, the corresponding guided drill is used. The tip of this drill, which has been specifically designed for this purpose, follows the pilot bore, which it widens until the planned diameter has been reached. The blunt tip automatically limits the depth of the cavity.

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During the next step, the depth bur extends the bore to the planned depth and width (the special design of the depth bur limits the depth of the bore here too). In order to achieve additional safety, depth markings are placed 2 mm apart.

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When conical implants are used, the final boring simply needs to be carried out with the corresponding drill of your implant system.

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- Ideally suited to prepare the jaw bone for the insertion of all commonly used implants
- Due to the special design, the guided and tri-spade drills follow the pre-defined pilot bore
- Use of only one drill sleeve per bore → cost effectiveness, cost reduction
- The innovative design ensures easy handling and very safe use
- Can be used for all three-dimensional X-ray images
- Additional optical control due to depth markings

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General instructions:

- All products delivered are unsterile therefore, before initial and each further potential application to be treated (cleaning / disinfection / sterilisation).

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