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- Ridge widening of the mandible by selective spreading and angle modulation of the implant bed
- Controlled, minimally invasive procedures
- Maximal conservation of existing bone substance
- Suitable for preliminary implant cavity preparation for all current implant systems
- Optimal in combination with Split-Control

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General instructions:

- All products delivered are unsterile therefore, before initial and each further potential application to be treated (cleaning / disinfection / sterilisation).

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Art.-No.: BHO00

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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¹ 204=RA, 206=RA XL, 314=FG

² Largest working part diameter in 1/10 mm

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Horizontal-Control is an Axial Angle Adjustment System that, with its completely new and unique instruments, allows simple and minimally invasive ridge widening. This is achieved because of an angle modulation of the implant axis given by the actual bone supply by erecting and displacing of the cortical plate, at the same time condensing the bone structure. In this way, the mandible can be optimally prepared for the use of all current implant systems. An important advantage of this innovative method is the non-traumatic procedure with special angled hand instruments. For special cases, Horizontal-Control is ideal when combined with the bone spreading instruments from the Split-Control system.

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An indication for the use of Horizontal-Control is when "pear-shaped" reabsorbed or traumatized sections of the lower alveolar ridge exist where conventional implantation without preliminary Bone Management® would not be possible.

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
First, the laser marked twist drill produces a sub-dimensional pilot hole which determines the depth and direction of the implant to be inserted (take notice of diameter and depth of the drilling)

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In the next step, the tapered expansion burs prepare a bone cavity appropriate for the subsequently used angle modulators.

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
The controlled used and implementation of a series of tapered angle modulators with increasing diameters, widens the implant cavity in the spongiosa and carefully condenses the neighboring bone tissue (Use the the corresponding angled attachment by screwing it onto the handle).

Successive alteration of the instrument alignment for angle modulation of the implant axis in the direction of the bone deficiency can erect detached segments, irrespective of the expansion of the alveolar ridge of up to 8 mm.

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
Should the existing jaw anatomy not allow enough expansion in the area of the implant cavity, a vertical osteotomy can be carried out (distal and mesial) internally without reflecting the periosteum. This causes a predetermined breaking point that allows the bone to break open in the following expansion

process. The bone casing is slightly displaced outwards so that the alveolar ridge can be widened to the necessary width.

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The subsequently inserted implant holds the expanded cavity open and stabilizes the cortical walls.

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In special indications when the alveolar ridge has to be spread excessively, a combined application with the spreading instruments from the Split-Control kit is beneficial.

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In this case, before the angle modulation, the deployment of a series of spreaders with increasing diameters causes the bone to be additionally spread and to be condensed.

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